

WWW.TEMPRANA.ORG



Options 1)

An initial three day home based parent training session with an approved Temprana Reflex Therapist.

The objective is to teach parents the specific application of a Temprana program for the child.

The application has to be applied 60 to 90 minutes per day and at least 4-6 days per week.

The parent training session is followed by a three hour consulting session by the Temprana Reflex Therapist every three months to monitor the child's development. The treatment plan will be adjusted accordingly.

Options 2)

Weekly 90 minute treatments with an approved Temprana Reflex Therapist.

The objective is that the Temprana Reflex Therapist carries out the child's individual and specific Temprana program on a weekly basis.

To devise an individual and specific Temprana program, a pre-consulting of 3 hours with a Temprana Reflex Therapist is necessary.

The program for children with special needs provides continual support to parents or other family members during the treatment period.

The objective is to provide an individual and specific service to each child to ensure the best level and quality in all sense of life is attained.

Temprana Reflex Therapy can be combined with all other methods of therapy and it is a method without any side effects.



Tel: 0034-933078972
sorensensistem@post.tele.dk
www.temprana.org
www.reflexologiafacial.es

WWW.TEMPRANA.ORG

KIDS TEMPRANA REFLEX
THERAPY

By Lone Sorensen Sorensensistem™

What is Temprana Reflex Therapy?

It is a natural and none invasive method for re-educating for brain functions for:

- Autism
- ADHD
- Hyperactivity
- Hypersensitivity
- Learning and behavior disabilities

International Institute
of Reflex Therapy



Temprana Reflex Therapy is a method to help re-educate the brain, Central Nerve System (CNS) and Hormonal system by manually stimulating areas and points related to the brain and CNS placed in the face, hands and feet. This is accomplished by using tools such as:

- Face Reflex therapy*
- Hand Reflex therapy*
- Foot Reflex therapy*

The basis for a specific program of stimulation for a child's individual situation is to analyze brain dysfunctions for:

- *Sensitivity*
- *Motor function*
- *Hyperactivity*
- *IQ*
- *Behavior*
- *Visual – and Auditory function.*
- *Learning, memory,*
- *concentration disabilities*
- *Immune system*

Including developing skills such as:

- *Hypo / Hyper sensitivity*
- *Social integration*
- *Language*
- *Communication*
- *Concentration*
- *Memory*
- *Learning*
- *Behavior*
- *Emotional difficulties*
- *Fine and gross motor function*
- *Co-ordination*

The success of Temprana Reflex Therapy is due to an individual analysis of dysfunctions in each child, the development of a tailored program of stimulation, and the carrying out of an intensive daily program.



Temprana means early intervention and the method is based on the most effective reflex therapy method worldwide.

Temprana Reflex Therapy should begin as early as possible. It is most effective when applied 8-10 hours weekly as team work between the family and a professional Temprana Reflex Therapist.

The International Institute of Temprana Reflex Therapy provides services for families with children with special needs including:

- Individual parent consulting*
- Home based parent courses*
- Individual therapy by a professional therapist*

The following modules have proven to be an effective recommendation for the application of Temprana Reflex Therapy sorenensistem TM.

